FORA 6

β-Ketone Test Strip

312-0099100-124 ver 1.0 2024/03

Warnings

- For *in vitro* diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this systems should handle everything that has comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your FORA 6 Duo / Plus / Connect / GTel Multi-functional meter Owner's Manual before you use this test strip.
- Use only FORA 6 β-Ketone test strips with FORA 6 Duo / Plus / Connect / GTel Multi-functional meter to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low bloodpressure or those who are in shock.
- Please do not use FORA 6 Duo / Plus / Connect / GTel Multi-functional monitoring system on critically ill patients. The collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological β-Ketone level might not be a true reflection. The following circumstances may apply: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemic, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA 6 β-Ketone test strips, when used together with FORA 6 Duo / Plus / Connect / GTel Multi-functional meter, allow your β-Ketone levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use ONLY heparin for anticoagulation of whole blood samples. Please do NOT use EDTA for anticoagulation.

Contents

- β-Ketone Test Strip
- Code strip
- Instruction for use

Limitations

- Hematocrit: The hematocrit level is limited to between 10% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: **This test strip must not be used for the testing of newborns.**
- This test strip is used for testing fresh capillary and venous.
- Altitude Effects: Altitudes up to 3,275 m (10,742 feet) do not affect test results.
- Please see Appendix: Summary of substances and concentrations in excess of limitation with interference.


Storage and Handling

IMPORTANT: Do not use the test strips if they have expired.

- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- Use each test strip immediately after taking it out of the vial or the individual foil packet. Close the vial immediately after taking out a strip. Keep the vial closed at all times.
- Store the test strips in their original vial ONLY. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Do not bend, cut, or alter the test strip.
- Store the test strips in a cool, dry place between 2 °C and 30 °C (35.6 °F and 86 °F) and 10% to 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.

Testing Your β-Ketone

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

The β-Ketone readings deliver plasma equivalent results and are displayed in millimoles of β-Ketone per liter of blood (mmol/L). The β-Ketone test measures β-Hydroxybutyrate (β-OHB), the most important of the three β-Ketone bodies in the blood. Normally, levels of β-OHB are expected to be less than 0.6 mmol/L. β-OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill. If your β-Ketone result is 0.0 mmol/L, repeat the β-Ketone test with new test strips. If the same message appears again or the result does not reflect how you feel, contact your healthcare professional. Follow your healthcare professional's advice before you make any changes to your diabetes medication programme. If your β-Ketone result is between 0.6 and 1.5 mmol/L, this may indicate development of a problem that could require medical assistance. Follow your healthcare professional's instructions. If your β-Ketone result is higher than 1.5 mmol/L, contact your healthcare professional promptly for assistance. You may be at risk of developing diabetic ketoacidosis (DKA).^{1,6}

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Upozornenie: Unusually high or low β-Ketone levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

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Table A summarized the substances tested at recommended concentrations which bias is less than 10% and demonstrated no interference with FORA Test Strip.

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Reference / Reference / Referencia

- [1] Schade DS, Eaton RP. Metabolic and clinical significance of ketosis. Special Topics in Endocrinology and Metabolism 1982;4:1-27.
- [2] Wiggam MI, O’Kane MJ, Harper R, Atkinson AB, Hadden DR, Trimble ER, Bell PM. Treatment of diabetic ketoacidosis using normalization of blood β-Hydroxybutyrate concentration as the endpoint of emergency management. Diabetes Care 1997;20:1347-52.
- [3] Harano Y, Kosugi K, Hyosu T, Suzuki M, Hidaka H, Kashiwagi A, Uno S, Shigeta Y. Ketone bodies as markers for Type 1 (insulin-dependent) diabetes and their value in the monitoring of diabetes control. Diabetologia 1984;26:343-8.
- [4] Ubukata E. Diurnal variation of blood β-Ketone bodies in insulin-dependent diabetes mellitus and noninsulindependent diabetes mellitus patients: The relationship to serum C-peptide immunoreactivity and free insulin. Ann Nutr Metab 1990;34:333-42.
- [5] Luzi L, Barrett EJ, Groop LC, Ferrannini E, DeFronzo RA. Metabolic effects of low-dose insulin therapy on glucose metabolism in diabetic ketoacidosis. Diabetes 1988;37:1470-77.
- [6] Hale PJ, Crase J, Natrass M. Metabolic effects of bicarbonate in the treatment of diabetic ketoacidosis. Br Med J 1984;289:1035-8.

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ForaCare Suisse AG
MedNet EC-REP GmbH